



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Wheaton sample vials, amber glass, without
closure, 8 ml PU 200 pcs**

**€184.84
plus VAT &
Shipping**

Product Images



Description

High-quality sample vials made of borosilicate glass resistant to temperature changes

Autoclavable

Packed in trays of 100 pieces each, dust-free.

Amber glass protects light-sensitive samples and maintains purity of contents.

Supplied without screw caps.

Additional Information

No.	WD-0710
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0161
EAN	4058072427061
Transport temperature	Room temperature
Volume	8 ml
Color bottle	brown
Height	60 mm
Floor outside	flat
Floor inside	flat
Thread	15-425 (US)
Neck size	far
DM bottle	17 mm
DM Closure	15 mm
Material bottle	Borosilicate glass (Wheaton 320. USP standard type 1)
TBST MAX bottle	510 C°

