



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Shorty-Vials, clear glass, 4 ml without closure, 200 pcs

€152.22
plus VAT &
Shipping

Product Images



Description

Low vials made of clear glass, without screw cap.

The wider shape results in greater stability and low dead volume above the sample material.

Temperature-resistant borosilicate glass Wheaton 33.

Packed dust-free in a practical Lab-File box.

Additional Information

| | |
|------------------------|--|
| No. | WD-0712 |
| Manufacturer (Brand) | Wheaton |
| old neoLab article no. | 9-0164 |
| EAN | 4058072427085 |
| Transport temperature | Room temperature |
| Volume | 4 ml |
| Color bottle | transparent (clear) |
| Height | 38 mm |
| Floor outside | flat |
| Floor inside | flat |
| DM bottle | 17 mm |
| DM Closure | 15 mm |
| Thread | 15-425 (US) |
| Neck size | far |
| Material bottle | Borosilicate glass (Wheaton 33. ASTM standard type 1 class A. USP standard type 1) |
| TBST MAX bottle | 505 C° |

