



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Shorty-Vials, clear glass, without closure, 6 ml PU 200 pcs. No.9-0166

€160.38
plus VAT &
Shipping

Product Images



Description

Low vials made of clear glass, without screw cap.

The wider shape results in greater stability and low dead volume above the sample material.

Temperature-resistant borosilicate glass Wheaton 33.

Packed dust-free in a practical Lab-File box.

Additional Information

No.	WD-0713
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0165
EAN	4058072427092
Transport temperature	Room temperature
Volume	6 ml
Color bottle	transparent (clear)
Height	40 mm
Floor outside	flat
Floor inside	flat
Thread	15-425 (US)
Neck size	far
DM bottle	19 mm
DM Closure	15 mm
Material bottle	Borosilicate glass (Wheaton 33. ASTM standard type 1 class A. USP standard type 1)
TBST MAX bottle	505 C°

