



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton sample vial 1.8 ml brown, with screw cap, 200 pcs./pack

**€256.88
plus VAT &
Shipping**

Product Images



Description

Brown Wheaton borosilicate glass sample vials 33 including autoclavable screw cap with rubber/Teflon® septum. Comply with ASTM Standard Class A/USP Class I.

Additional Information

No.	WD-0695
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0117
EAN	4058072426958
Transport temperature	Room temperature
Volume	1.8 ml
Color bottle	brown
Height	35 mm
Floor outside	flat
Floor inside	flat
Thread	08-425 (US)
Quantity components	two
Neck size	narrow
Closure / lid properties	closed
DM bottle	12 mm
DM Neck	8 mm
Color closure/ lid	black
Material bottle	Borosilicate glass (Wheaton 320. USP standard type 1)
Material closure/ lid	Polypropylene (PP)

