



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Wheaton sample vial brown 8 ml, cap+seal, 144
pcs./pack**

**€158.54
plus VAT &
Shipping**

Product Images



Description

Brown Wheaton 320 borosilicate glass sample vials for light sensitive sample material. The glass meets ASTM Type I Checkout B requirements.

Sample vials and caps can be sterilized and stored at low temperatures.

Screw caps with white rubber inserts.

Volume: 8.0 ml

Additional Information

No.	WD-0730
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0180
EAN	4058072427139
Transport temperature	Room temperature
Volume	8 ml
Color bottle	brown
Height	63 mm
Floor outside	flat
Floor inside	flat
Thread	15-425 (US)
Quantity components	a
Neck size	narrow
Closure / lid properties	closed
DM bottle	17 mm
DM Neck	15 mm
Color closure/ lid	black
Material bottle	Borosilicate glass (Wheaton 320. USP standard type 1)
Material closure/ lid	Polypropylene (PP)

