



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® short thread vials ND9, 1.5 ml clear
glass, silanized, wide opening, 100 pcs./pack**

**€44.60
plus VAT &
Shipping**



Description

Due to their technical geometry, the bottles can be used on all common autosamplers. The wide opening allows easy filling with viscous substances.

Additional Information

No.	7-0674
Manufacturer (Brand)	neochrom
EAN	4058072077099
Transport temperature	Room temperature
Volume	1.5 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat
Floor inside	flat
DM bottle	11.6 mm
DM Neck	9 mm
DM Closure	9 mm
Neck size	far
Material bottle	Borosilicate glass
ND (nominal diameter)	ND9

