



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® 2 ml screw vials, clear glass, with  
title block, thread 9-425, blue cap, red  
PTFE/white silicone; Packing unit: 100 pcs**

**€49.50  
plus VAT &  
Shipping**

## Product Images



## Description

The neochrom Combi Kits consisting of 2 ml vials with 9-425 thread and the corresponding lids. All combination kits are suitable for most applications with Agilent autosamplers. Kits are unmounted and consist of 100 vials, plus 100 matching caps. \* slotted caps recommended for UHPLC

## Additional Information

No.	EC-1191
Manufacturer (Brand)	neochrom
EAN	4058072080280
Transport temperature	Room temperature
Volume	2 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat
Floor inside	flat
DM Hole	6 mm
Thread	09-425 (US)
Quantity components	two
Neck size	narrow
Closure / lid properties	pierceable
DM bottle	12 mm
DM Neck	8 mm
Color closure/ lid	blue
DM Closure	10 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)

