



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® HPLC/LC-MS cert. 2 ml screw vial,
clear glass, with title block, thread 9-425, blue
cap, red PTFE/white silicone with slot*; Packing
unit: 100 pcs**

**€74.80
plus VAT &
Shipping**

Product Images



Description

The neochrom Combi Kits consisting of 2 ml vials with thread 9-425 and the corresponding lids. All combi kits are certified for HPLC/LC-MS applications for Waters autosamplers. Kits are unmounted and consist of 100 vials, plus 100 matching caps. * slotted caps recommended for Waters UPLC™ and UHPLC

Additional Information

No.	EC-1186
Manufacturer (Brand)	neochrom
EAN	4058072080150
Transport temperature	Room temperature
Volume	2 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat
Floor inside	flat
DM Hole	6 mm
Thread	09-425 (US)
Quantity components	two
Neck size	narrow
Closure / lid properties	pierceable
DM bottle	12 mm
DM Neck	8 mm
Color closure/ lid	blue
DM Closure	10 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)

