



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® 2 ml screw vials, clear glass, with
title block, mounted with 9 mm blue cap, red
PTFE/white silicone; Packing unit: 100 pcs**

**Price on
request
plus VAT &
Shipping**



Description

The neochrom Combi Kits consist of 100 x 2 ml vials with 9-425 thread and the corresponding pre-mounted lids. All combi kits are suitable for most applications with common autosampler models.

Additional Information

No.	EC-1205
Manufacturer (Brand)	neochrom
EAN	4058072080426
Transport temperature	Room temperature
Volume	2 ml
Color bottle	transparent (clear)
Height	32 mm
Quantity components	two
Floor outside	flat
Floor inside	flat
DM bottle	12 mm
DM Neck	8 mm
DM Hole	6 mm
DM Closure	10 mm
Color closure/ lid	blue
Thread	09-425 (US)
Neck size	narrow
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	pierceable

