



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® HPLC/LC-MS cert. 2 ml screw vial,
amber glass, with title block, thread 9-425, blue
cap, red PTFE/white silicone; Packing unit: 100
pcs**

**€64.00
plus VAT &
Shipping**

Product Images



Description

The neochrom Combi Kits consist of 2 ml vials with thread 9-425 and the corresponding lids. All combi kits are certified for HPLC/LC-MS applications for Agilent autosamplers. Kits are unmounted and consist of 100 vials, plus 100 matching caps. * slotted caps recommended for UHPLC

Additional Information

No.	EC-1181
Manufacturer (Brand)	neochrom
EAN	4058072078904
Transport temperature	Room temperature
Volume	2 ml
Color bottle	brown
Height	32 mm
Quantity components	two
Floor outside	flat
Floor inside	flat
DM bottle	12 mm
DM Neck	8 mm
DM Hole	6 mm
DM Closure	10 mm
Color closure/ lid	blue
Thread	09-425 (US)
Neck size	narrow
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	pierceable

