



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® snap ring micro vial ND9 0.3 ml, TPX highly transparent, 100 pcs./pack

Special Price

€21.00 was

€31.50

**plus VAT &
Shipping**



Description

Snap ring bottles ND11

Made of TPX with integrated glass micro insert.
TPX is a plastic with crystal clear transparency and good chemical resistance.

Dimensions: 32 x 11.6 mm Ø

Additional Information

No.	7-0722
Manufacturer (Brand)	neochrom
EAN	4058072029777
Transport temperature	Room temperature
Volume	0.25 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	conical. pointed . with stand ring
Floor inside	conical. pointed
Neck size	far
DM bottle	11.6 mm
DM Neck	11 mm
DM Closure	11 mm
Material bottle	Polymethylpentene (PMP. TPX®)
for the autosampler model	Agilent 1050 Agilent 1090 Agilent 1100 Agilent 1200 Agilent 7693A Agilent 79855(A) Agilent CTC HTS/ HTC PAL Alcott 719 AL Alcott 719 D/ D-PCS Alcott 738 Antec-Leyden AS 100 Antek 736 Unisampler Antek 738 Beckmann 501 Beckmann 502/ 502e Beckmann 507/ 507e

