



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® Snap ring vials PP 0.3 ml ND11,  
brown, 32 x 11.6 mm, 100 pcs./pack**

**€22.50  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Snap ring bottles ND11

Made of PP with integrated micro insert

Dimensions: 32 x 11.6 mm Ø

## Additional Information

---

No.	7-0724
Manufacturer (Brand)	neochrom
EAN	4058072029784
Transport temperature	Room temperature
Volume	0.25 ml
Color bottle	brown
Height	32 mm
Floor outside	conical. pointed . with stand ring
Floor inside	conical. pointed
Neck size	far
DM bottle	11.6 mm
DM Neck	11 mm
DM Closure	11 mm
Material bottle	Polypropylene (PP)
for the autosampler model	Agilent 1050 Agilent 1090 Agilent 1100 Agilent 1200 Agilent 7693A Agilent 79855(A) Agilent CTC HTS/ HTC PAL Alcott 719 AL Alcott 719 D/ D-PCS Alcott 738 Antec-Leyden AS 100 Antek 736 Unisampler Antek 738 Beckmann 501 Beckmann 502/ 502e Beckmann 507/ 507e

