



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neochrom® TopSert snap ring vial ND11, 0.2 ml, TPX, with glass micro insert, 100 pcs./pack

Special Price

**€45.00** was

**€67.50**

**plus VAT &  
Shipping**



## Description

---

Snap ring bottles ND11

Made of TPX with integrated glass micro insert.  
TPX is a plastic with crystal clear transparency and good chemical resistance.

Dimensions: 32 x 11,6 mm Ø

## Additional Information

---

|                           |  |
|---------------------------|--|
| No.                       | 7-0627   |
| Manufacturer (Brand)      | neochrom   |
| EAN                       | 4058072028985  |
| Transport temperature     | Room temperature   |
| Volume                    | 0.2 ml   |
| Color bottle              | transparent (clear)  |
| Height                    | 32 mm  |
| Floor outside             | conical. pointed . with stand ring   |
| Floor inside              | conical. pointed   |
| Neck size                 | far  |
| DM bottle                 | 11.6 mm  |
| DM Neck                   | 11 mm  |
| DM Closure                | 11 mm  |
| Material bottle           | Polymethylpentene (PMP. TPX®)  |
| for the autosampler model | Agilent 1050 Agilent 1090 Agilent 1100 Agilent 1200 Agilent 7693A Agilent 79855(A) Agilent CTC HTS/ HTC PAL Alcott 719 AL Alcott 719 D/ D-PCS Alcott 738 Antec-Leyden AS 100 Antek 736 Unisampler Antek 738 Beckmann 501 Beckmann 502/ 502e Beckmann 507/ 507e |

