



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® snap ring vials 1.5 ml ND11, wide opening, clear glass, 100 pcs./pack

**€12.80
plus VAT &
Shipping**



Description

Snap ring bottles ND11

Glass 1st hydrolytic class.

The glass bottles can also be closed with a standard 11 mm flanged cap, as the two snap ring lips together have the same height as a rolled rim.

Dimensions: 32 x 11.6 mm Ø

Additional Information

No.	7-0748
Manufacturer (Brand)	neochrom
EAN	4058072029999
Transport temperature	Room temperature
Volume	1.5 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat
Floor inside	flat
Neck size	far
DM bottle	11.6 mm
DM Neck	11 mm
DM Closure	11 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent 1050 Agilent 1090 Agilent 1100 Agilent 1200 Agilent 7693A Agilent 79855(A) Agilent CTC HTS/ HTC PAL Alcott 719 AL Alcott 719 D/ D-PCS Alcott 738 Antec-Leyden AS 100 Antek 736 Unisampler Antek 738 Beckmann 501 Beckmann 502/ 502e Beckmann 507/ 507e

