



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® snap ring vial amber glass, 2 ml,
opening 11 mm, 12 x 32 mm, with title block**

**€14.80
plus VAT &
Shipping**



Description

The neochrom snap ring vials ND11, vol. 2 ml, 12 x 32 mm are made of borosilicate glass of the 1st hydrolytic class.

The cooling point is 565 °C.

Suitable for AGILENT™ (HP) and Perkin Elmer autosamplers.

Please order matching snap ring caps separately.

Additional Information

No.	EC-1054
Manufacturer (Brand)	neochrom
EAN	4058072080044
Transport temperature	Room temperature
Volume	2 ml
Color bottle	brown
Height	32 mm
Floor outside	flat
Floor inside	flat
DM bottle	12 mm
DM Neck	11 mm
DM Closure	11 mm
Neck size	far
Material bottle	Borosilicate glass (1st hydrolytic class)

