



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® EPA thread vials ND24, 20.0 ml clear
glass, 57 x 27.5 mm, 100 pcs./pack**

**€41.20
plus VAT &
Shipping**



Description

Threaded bottles (EPA) made of glass of 1st hydrolytic class with thread 24-400

Pack unit: 100 pcs/pack.

Additional Information

No.	7-0851
Manufacturer (Brand)	neochrom
EAN	4058072037192
Transport temperature	Room temperature
Volume	20 ml
Color bottle	transparent (clear)
Height	57 mm
Floor outside	flat
Floor inside	flat
DM bottle	27.5 mm
DM Closure	24 mm
Thread	24-400 (US)
Neck size	far
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent 7695A Agilent AQUATek 70 Agilent Archon Pluge + Trap Agilent Tekmar SOLATek72 Analytik Jena multi N/C 3000 (TOC) Dimatec Dimatoc 200 Dimatec Dimatoc 300 Dimatec Dimatoc 400 Dimatec Dimatoc 2000 Dionex ASE 200 EST Analytical Archon P/T EST Analytic

