



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® EPA Thread Vials ND24, 40,0 ml  
clear glass, 95 x 27,5 mm, Thread 24-400 1st  
hydrolytic class ( Pack á 100 pieces )**

**€43.60  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

**Threaded bottles (EPA) made of glass of 1st hydrolytic class with thread 24-400**

Pack unit: 100 pcs/pack.

## Additional Information

---

No.	7-0855
Manufacturer (Brand)	neochrom
EAN	4058072037239
Transport temperature	Room temperature
Volume	40 ml
Color bottle	transparent (clear)
Height	95 mm
Floor outside	flat
Floor inside	flat
Thread	24-400 (US)
Neck size	far
DM bottle	27.5 mm
DM Closure	24 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent 7695A Agilent AQUATek 70 Agilent Archon Pluge + Trap Agilent Tekmar SOLATek72 Analytik Jena multi N/C 3000 (TOC) Dimatec Dimatoc 200 Dimatec Dimatoc 300 Dimatec Dimatoc 400 Dimatec Dimatoc 2000 Dionex ASE 200 EST Analytical Archon P/T EST Analytic

