



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® fine thread vials ND18, amber glass,
10 ml, 100 pcs./pack**

**€48.40
plus VAT &
Shipping**



Description

Fine thread glassbottles of 1st hydrolytic class, with rounded bottom.

Additional Information

No.	7-0841
Manufacturer (Brand)	neochrom
EAN	4058072030759
Transport temperature	Room temperature
Volume	8 ml
Color bottle	brown
Height	46 mm
Floor outside	rounded transition to the bottle bottom (HS bottom)
Floor inside	rounded transition to the bottle bottom (HS bottom)
DM bottle	22.5 mm
DM Closure	18 mm
Thread	Fine thread
Neck size	far
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent CTC Combi PAL AI Headspace Carlo Erba HS800 Carlo Erba HS850 CTC (LEAP) Combi PAL/ GC PAL. 32 Pos. Tray CTC (LEAP) Combi PAL SMPE Mode. 32 Pos. Tray CTC (LEAP) HTX PAL/ HTC PAL/ HTS PAL. 32 Pos. Tray CTC Combi-xt Headspace Option CTC Combi-xt Liq

