



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® fine thread vials ND18, clear glass,
20 ml, 100 pcs./pack**

**€34.20
plus VAT &
Shipping**

Product Images



Description

Fine thread glassbottles of 1st hydrolytic class, with rounded bottom.

Additional Information

No.	7-0842
Manufacturer (Brand)	neochrom
EAN	4058072037116
Transport temperature	Room temperature
Volume	18 ml
Color bottle	transparent (clear)
Height	75.5 mm
Floor outside	rounded transition to the bottle bottom (HS bottom)
Floor inside	rounded transition to the bottle bottom (HS bottom)
Thread	Fine thread
Neck size	far
DM bottle	22.5 mm
DM Closure	18 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent CTC Combi PAL AI Headspace Carlo Erba HS800 Carlo Erba HS850 CTC (LEAP) Combi PAL/ GC PAL. 32 Pos. Tray CTC (LEAP) Combi PAL SMPE Mode. 32 Pos. Tray CTC (LEAP) HTX PAL/ HTC PAL/ HTS PAL. 32 Pos. Tray CTC Combi-xt Headspace Option CTC Combi-xt Liq

