



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® fine thread vials ND18, amber glass,  
20 ml, 100 pcs./pack**

**€52.20  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Fine thread glassbottles of 1st hydrolytic class, with rounded bottom.

## Additional Information

---

No.	7-0843
Manufacturer (Brand)	neochrom
EAN	4058072037123
Transport temperature	Room temperature
Volume	18 ml
Color bottle	brown
Height	75.5 mm
Floor outside	rounded transition to the bottle bottom (HS bottom)
Floor inside	rounded transition to the bottle bottom (HS bottom)
DM bottle	22.5 mm
DM Closure	18 mm
Thread	Fine thread
Neck size	far
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent CTC Combi PAL AI Headspace Carlo Erba HS800 Carlo Erba HS850 CTC (LEAP) Combi PAL/ GC PAL. 32 Pos. Tray CTC (LEAP) Combi PAL SMPE Mode. 32 Pos. Tray CTC (LEAP) HTX PAL/ HTC PAL/ HTS PAL. 32 Pos. Tray CTC Combi-xt Headspace Option CTC Combi-xt Liq

