



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® Autosampler Vials 1.2 ml, clear glass, 32 x 12 mm ND8, 100 pcs./pack

€30.70
plus VAT &
Shipping

Product Images



Description

Autosampler vials 1.2 ml with 8-425 thread,
made of borosilicate glass type 1.
height 32 mm, Ø 12 mm.

Additional Information

No.	2-7200
Manufacturer (Brand)	neochrom
EAN	4058072008987
Transport temperature	Room temperature
Volume	1.2 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat
Floor inside	flat
Thread	08-425 (US)
Neck size	narrow
DM bottle	12 mm
DM Neck	8 mm
DM Closure	8 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	KA

