



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® Autosampler Vials 1.2 ml, amber glass, 32 x 12 mm ND8, 100 pcs./pack

€32.60
plus VAT &
Shipping

Product Images



Description

Autosampler vials 1.2 ml with 8-425 thread,
made of borosilicate glass type 1.
height 32 mm, Ø 12 mm.

Additional Information

No.	2-7201
Manufacturer (Brand)	neochrom
EAN	4058072008994
Transport temperature	Room temperature
Volume	1.2 ml
Color bottle	brown
Height	32 mm
Floor outside	flat
Floor inside	flat
DM bottle	12 mm
DM Neck	8 mm
DM Closure	8 mm
Thread	08-425 (US)
Neck size	narrow
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	KA

