



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® thread vials ND8, 1.1 ml clear glass
32 x 11.6 mm, conical, 100 pcs./pack**

**€52.20
plus VAT &
Shipping**



Description

Standard sample vial with 8-425 thread for GC and HPLC.

For almost all common autosamplers.

Made of glass of the 1st hydrolytic class.

The narrow opening requires micro-inserts with a diameter of 5 mm.

Additional Information

No.	7-0659
Manufacturer (Brand)	neochrom
EAN	4058072029234
Transport temperature	Room temperature
Volume	1 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	conical. pointed
Floor inside	conical. pointed
DM bottle	11.6 mm
DM Neck	8 mm
DM Closure	8 mm
Thread	08-425 (US)
Neck size	narrow
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Al 42 vial tray Beckmann 507/ 507e Beckmann Marathon Beckmann Promis Beckmann Triathlon Standard Tray Carlo Erba A200LC Carlo Erba AS100 Carlo Erba AS200S Carlo Erba AS300 Cecil Instruments AutoQuest CTC (LEAP) HTX PAL/ HTC PAL/ HTS PAL. 54/98 Pos. Tray C

