



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® wide-neck threaded vials amber
glass with title block, 2 ml, thread 10 mm, 12 x 32
mm**

**€16.00
plus VAT &
Shipping**

Product Images



Description

The neochrom threaded vials ND10, vol. 2 ml, 12 x 32 mm are made of borosilicate glass of the 1st hydrolytic class. The cooling point is 565 °C.

The vials are suitable for Bruker (formerly Varian) autosamplers.

Please order matching screw caps separately.

Additional Information

No.	EC-1111
Manufacturer (Brand)	neochrom
EAN	4058072080228
Transport temperature	Room temperature
Volume	2 ml
Color bottle	brown
Height	32 mm
Floor outside	flat
Floor inside	flat
Thread	10-425 (US)
Neck size	far
DM bottle	12 mm
DM Closure	10 mm
Material bottle	Borosilicate glass (1st hydrolytic class)

