



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® thread vials ND10, clear glass, 1.5 ml, 100 pcs/pack

**€13.50
plus VAT &
Shipping**



Description

Standard sample vials made of clear or amber glass with wide opening.

1.5 ml capacity, thread 10-425.

Made of 1st hydrolytic class glass.

Dimensions: 32 x 11.6 mm Ø.

Additional Information

No.	EC-1308
Manufacturer (Brand)	neochrom
EAN	4058072076894
Transport temperature	Room temperature
Volume	2 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat
Floor inside	flat
DM bottle	12 mm
DM Closure	10 mm
Thread	10-425 (US)
Neck size	far
Material bottle	Borosilicate glass (1st hydrolytic class)

