



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® thread vials ND13, 4,0 ml clear glass
45 x 14,7 mm, thread 13-425, text field + filling
marks 100 pcs./pack**

**€17.80
plus VAT &
Shipping**

Product Images



Description

Standard sample vials made of clear or amber glass with wide opening.

4 ml capacity, thread 13-425.

Made of 1st hydrolytic class glass.

Dimensions: 45 x 14.7 mm Ø.

Additional Information

No.	7-0765
Manufacturer (Brand)	neochrom
EAN	4058072030148
Transport temperature	Room temperature
Volume	4.1 ml
Color bottle	transparent (clear)
Height	45 mm
Floor outside	flat
Floor inside	flat
Thread	13-425 (US)
Neck size	far
DM bottle	14.7 mm
DM Closure	13 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent G1888A Agilent HS 7694 Agilent 6850. 22 pos. tray Agilent 7693A Alcott 719 AL Beckmann Triathlon LSV tray Bruker LC51 Cambridge Scientific Instruments 300 Series Cambridge Scientific Instruments 205 Series Cecil Instruments AutoQuest Dionex AS 40

