



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® thread vials ND13, 4.0 ml amber
glass with writing field and filling markings, 100
pcs./pack**

**€18.90
plus VAT &
Shipping**



Description

Standard sample vials made of clear or amber glass with wide opening.

4 ml capacity, thread 13-425.

Made of 1st hydrolytic class glass.

Dimensions: 45 x 14.7 mm Ø.

Additional Information

No.	7-0767
Manufacturer (Brand)	neochrom
EAN	4058072030162
Transport temperature	Room temperature
Volume	4.1 ml
Color bottle	brown
Height	45 mm
Floor outside	flat
Floor inside	flat
DM bottle	14.7 mm
DM Closure	13 mm
Thread	13-425 (US)
Neck size	far
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent G1888A Agilent HS 7694 Agilent 6850. 22 pos. tray Agilent 7693A Alcott 719 AL Beckmann Triathlon LSV tray Bruker LC51 Cambridge Scientific Instruments 300 Series Cambridge Scientific Instruments 205 Series Cecil Instruments AutoQuest Dionex AS 40

