



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® TopSert TPX short thread vial, 32 x 11,6mm, with integrated 0,2ml glass micro insert, 15mm tip, brown (Pack á 100 pieces)

**€68.50
plus VAT &
Shipping**

Product Images



Description

Main applications for plastic bottles are atomic absorption, water and protein analysis, capillary electrophoresis, ion chromatography, etc.

- TPX has crystal clear transparency and good chemical resistance

Dimensions: 32 x 11.6 mm

Additional Information

No.	7-0682
Manufacturer (Brand)	neochrom
EAN	4058072029456
Transport temperature	Room temperature
Volume	0.2 ml
Color bottle	brown
Height	32 mm
Floor outside	conical. pointed . with stand ring
Floor inside	conical. pointed
DM bottle	11.6 mm
DM Closure	9 mm
Thread	Short thread
Neck size	far
Material bottle	Polymethylpentene (PMP. TPX®)
for the autosampler model	Agilent 1050 Agilent 1090 Agilent 1100 Agilent 1200 Agilent 5880 Agilent 5890 Agilent 6850. 27 Pos. Tray Agilent 6890 Agilent 7673A Agilent 7683A Agilent 79855(A) Agilent CTC GC PAL Agilent CTC HTS/ HTC PAL AI 42 vial tray AI 60 vial tray AI CTC A200S AIM

