



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® short thread vials ND9 0.3 ml, PP
transparent, 100 pcs./pack**

**€17.30
plus VAT &
Shipping**



Description

Main applications for plastic vials are atomic absorption, water and protein analysis, capillary electrophoresis, ion chromatography, etc.

- 1.5 ml short thread vials with slightly concave bottom:
applicable in HPLC as an alternative to silanized glass vials for pH-sensitive and aqueous samples

Dimensions: 32 x 11.6 mm

Additional Information

No.	7-0684
Manufacturer (Brand)	neochrom
EAN	4058072029470
Transport temperature	Room temperature
Volume	0.25 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	conical. pointed . with stand ring
Floor inside	conical. pointed
Thread	Short thread
Neck size	far
DM bottle	11.6 mm
DM Closure	9 mm
Material bottle	Polypropylene (PP)
for the autosampler model	Agilent 1050 Agilent 1090 Agilent 1100 Agilent 1200 Agilent 5880 Agilent 5890 Agilent 6850. 27 Pos. Tray Agilent 6890 Agilent 7673A Agilent 7683A Agilent 79855(A) Agilent CTC GC PAL Agilent CTC HTS/ HTC PAL AI 42 vial tray AI 60 vial tray AI CTC A200S AIM

