



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® Micro Vials 0.3 ml, cylindrical w.
inner cone, made of PP, 100 pcs/pack**

**€14.50
plus VAT &
Shipping**

Product Images



Description

Microsample vessels for HPLC and GC autosamplers, suitable for most standard autosamplers.
Neck diameter 10.8 mm. Made of polypropylene, milky transparent, self-standing, autoclavable, suitable for freezing.

Additional Information

No.	2-1012
Manufacturer (Brand)	neochrom
EAN	4058072072216
Volume	0.3 ml
Color bottle	transparent (milky)
Height	31.8 mm
Locking micro insert	at the bottleneck
Floor outside	flat
Floor inside	conical. pointed
DM Neck	10.8 mm
Material bottle	Polypropylene (PP)
for the autosampler model	KA

