



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® ND20 headspace vials, 5 ml clear  
glass, 38.2 x 22 mm, round bottom, 100  
pcs./pack**

**€26.80  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

The instruments on the market require a wide variety of bottles.

To help you identify the correct bottle, we have provided instrument suitability information in our product overview.

All bottles are 1st hydrolytic class glass.

Headspace bottles differ in the following characteristics:

- Volume
- rounded or flat bottom
  - A rounded bottom is more stable at high internal pressure during tempering. In addition, the bottle slides more easily into the heating block when inserted there by means of magnets.
  - A flat bottom may be necessary if the bottles have to run slightly downwards in the instrument.
- beveled headspace rim
  - A beveled headspace rim may be necessary for some special closure systems (Perkin Elmer).
- Neck length (instrument-dependent)
- Clear / amber glass
- with / without title block and filling marks

## Additional Information

No.	7-0788
Manufacturer (Brand)	neochrom
EAN	4058072030346
Transport temperature	Room temperature
Volume	5 ml
Color bottle	transparent (clear)
Height	38.2 mm
Floor outside	rounded transition to the bottle bottom (HS bottom)
Floor inside	rounded transition to the bottle bottom (HS bottom)
DM bottle	22 mm
DM Closure	20 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
Type bottle rim	Headspace edge
for the autosampler model	Alcott 719 AL Beckmann Triathlon Standard Tray Knauer K-3800 (Basic Marathon) Knauer Smartline K-3950 Perkin Elmer HS6 Perkin Elmer HS40/ HS100/ HS101 Perkin Elmer ISS-100. 85 vial tray Perkin Elmer ISS-200. 85 vial tray Perkin Elmer ISS-225. 25 vial tray

