



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® Headspace vials, 20 ml, clear glass,  
beaded rim, 75.5 x 22.5 mm, rounded bottom,  
type Varian, 100 pcs./pack**

**€29.90  
plus VAT &  
Shipping**



## Description

---

The instruments on the market require a wide variety of bottles.

To help you identify the correct bottle, we have provided information on instrument suitability in our product overview. All bottles are made of 1st hydrolytic class glass.

Headspace bottles differ in the following characteristics:

- Volume
- rounded or flat bottom
  - A rounded bottom is more stable at high internal pressure during tempering. In addition, the bottle slides more easily into the heating block when inserted there by means of magnets.
  - A flat bottom may be necessary if the bottles have to run slightly downwards in the instrument.
- beveled headspace rim
  - A beveled headspace rim may be necessary for some special closure systems (Perkin Elmer).
- Neck length (instrument-dependent)
- Clear / amber glass
- with / without title block and filling marks

## Additional Information

No.	7-0792
Manufacturer (Brand)	neochrom
EAN	4058072030384
Transport temperature	Room temperature
Volume	20 ml
Color bottle	transparent (clear)
Height	75.5 mm
Floor outside	rounded transition to the bottle bottom (HS bottom)
Floor inside	rounded transition to the bottle bottom (HS bottom)
DM bottle	22.5 mm
DM Closure	20 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
Type bottle rim	DIN rolled edge
for the autosampler model	Agilent CTC Combi PAL AI Headspace Carlo Erba HS800 Carlo Erba HS850 CTC (LEAP) Combi PAL/ GC PAL. 32 Pos. Tray CTC (LEAP) Combi PAL SMPE Mode. 32 Pos. Tray CTC (LEAP) HTX PAL/ HTC PAL/ HTS PAL. 32 Pos. Tray CTC Combi-xt Headspace Option CTC Combi-xt Liq

