



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® rolled rim vials 0.8 ml clear glass,
30 x 8.2 mm, flat bottom, 100 pcs./pack**

**€12.80
plus VAT &
Shipping**

Product Images



Description

Sample vials with DIN rolled rim ND8 in various volumes and shapes.

For almost all common autosamplers.

Made of 1st hydrolytic class glass.

They can be closed with 8 mm aluminum crimp caps, with 9 mm PE caps or with 8 mm push-on caps.

Additional Information

No.	7-0600
Manufacturer (Brand)	neochrom
EAN	4058072028732
Transport temperature	Room temperature
Volume	0.8 ml
Color bottle	brown
Height	30 mm
Floor outside	flat
Floor inside	flat
DM bottle	8.2 mm
DM Closure	8 mm
Neck size	narrow
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	AI 60 vial tray Carlo Erba AS800. 60 vial tray Dynatech 60 vial tray Dynatech GC111 Dynatech GC311 Fisons AS800 60 vial tray LDC 713-60 Perkin Elmer AI-1 Perkin Elmer AS100/ AS100B Perkin Elmer AS2000/ AS2000B Perkin Elmer AS300 Perkin Elmer AS8300 Perkin

