




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® rolled rim vials 0.7 ml clear glass,
40 x 7 mm, flat bottom, 100 pcs./pack**

**€17.50
plus VAT &
Shipping**



Description

Sample vials with DIN rolled rim ND8 in various volumes and shapes.

For almost all common autosamplers.

Made of 1st hydrolytic class glass.

They can be closed with 8 mm aluminum crimp caps, with 9 mm PE caps or with 8 mm push-on caps.

Additional Information

No.	7-0601
Manufacturer (Brand)	neochrom
EAN	4058072028749
Transport temperature	Room temperature
Volume	0.8 ml
Color bottle	transparent (clear)
Height	40 mm
Floor outside	flat
Floor inside	flat
Neck size	narrow
DM bottle	7 mm
DM Closure	8 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent CTC A200S Agilent CTC GC PAL Agilent CTC HTS/ HTC PAL Beckmann 504 Beckmann Triathlon Micro-Tray Carlo Erba AS100 Carlo Erba AS200 Carlo Erba A200LC Carlo Erba AS200S Carlo Erba AS300 CTC (LEAP) Combi PAL/ GC PAL. 200 pos. tray CTC (LEAP) HTX PAL/

