



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® rolled rim vials 0.3 ml clear glass,  
31.5 x 5.5 mm, round bottom, 100 pcs./pack**

**€23.40  
plus VAT &  
Shipping**



## Description

---

Sample vials with DIN rolled rim ND8 in various volumes and shapes.

For almost all common autosamplers.

Made of 1st hydrolytic class glass.

They can be closed with 8 mm aluminum crimp caps, with 9 mm PE caps or with 8 mm push-on caps.

## Additional Information

---

No.	7-0603
Manufacturer (Brand)	neochrom
EAN	4058072028763
Transport temperature	Room temperature
Volume	0.3 ml
Color bottle	transparent (clear)
Height	31.5 mm
Floor outside	conical
Floor inside	around
DM bottle	5.5 mm
DM Closure	8 mm
Neck size	narrow
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent 1050. 34 pos. tray Agilent 1090. 34 pos. tray Agilent 7673A Agilent 7683A Agilent CTC HTS/ HTC PAL Agilent CTC GC PAL AI 60 vial tray Beckmann 501 Beckmann 502/ 502e Beckmann 507/ 507e Beckmann Triathlon LSV tray Carlo Erba AS100 Carlo Erba AS300

