



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® rolled rim vials 1.2 ml amber glass
40 x 8.2 mm, 1st hydr. class, 100 pcs./pack**

**€18.20
plus VAT &
Shipping**



Description

Sample vials with DIN rolled rim ND8 in various volumes and shapes.

For almost all common autosamplers.

Made of 1st hydrolytic class glass.

They can be closed with 8 mm aluminum crimp caps, with 9 mm PE caps or with 8 mm push-on caps.

Additional Information

No.	7-0606
Manufacturer (Brand)	neochrom
EAN	4058072028794
Transport temperature	Room temperature
Volume	1 ml
Color bottle	brown
Height	40 mm
Floor outside	flat
Floor inside	flat
Neck size	narrow
DM bottle	8.2 mm
DM Closure	8 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	AI 60 vial tray Carlo Erba AS800. 60 vial tray Dionex Gina 50 Dionex Summit ASI 100 micro tray 192 pos. Dionex UltiMate Nano/Cap/Micro WPS-3000 SL 216 (3x72) pos. rack 1.2 ml Dynatech 60 vial tray Dynatech LC2000 Fisons AS800 60 vial tray Gynkotek Gina 5

