



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657




neochrom® rolled rim micro vials 0.6 ml amber glass, conical, 40 x 7 mm, 100 pcs./pack

Special Price

€6.30 was

€9.45

**plus VAT &
Shipping**



Description

Sample vials with DIN rolled rim ND8 in various volumes and shapes.

For almost all common autosamplers.

Made of 1st hydrolytic class glass.

They can be closed with 8 mm aluminum crimp caps, with 9 mm PE caps or with 8 mm push-on caps.

Additional Information

No.	7-0608
Manufacturer (Brand)	neochrom
EAN	4058072028817
Transport temperature	Room temperature
Volume	0.6 ml
Color bottle	brown
Height	40 mm
Floor outside	conical. pointed
Floor inside	conical. pointed
Neck size	narrow
DM bottle	7 mm
DM Closure	8 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Agilent CTC A200S Agilent CTC GC PAL Agilent CTC HTS/ HTC PAL Beckmann 504 Beckmann Triathlon Micro-Tray Carlo Erba AS100 Carlo Erba AS200 Carlo Erba A200LC Carlo Erba AS200S Carlo Erba AS300 CTC A200S CTC A200LC Finnigan A200S Fisons AS100 Fisons AS200 F

