



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® rolled rim vials ND11, clear glass, 1.1 ml 32 x 11.6 mm, 1.hydr.class, 100 pcs./pack**

**€75.00  
plus VAT &  
Shipping**



## Description

---

Standard vials for GC and HPLC with rolled rim ND 11 for most autosamplers.

The wide opening allows easy filling.

Glass 1st hydrolytic class.

Vial dimensions: 32 x 11.6 mm (H x Ø)

## Additional Information

---

|                           |   |
|---------------------------|---|
| No.                       | 7-0625  |
| Manufacturer (Brand)      | neochrom  |
| EAN                       | 4058072028978   |
| Transport temperature     | Room temperature  |
| Volume                    | 1.5 ml  |
| Color bottle              | transparent (clear)   |
| Height                    | 32 mm   |
| Floor outside             | flat  |
| Floor inside              | conical. pointed  |
| DM bottle                 | 11.6 mm   |
| DM Closure                | 11 mm   |
| Neck size                 | far   |
| Material bottle           | Borosilicate glass (1st hydrolytic class)   |
| for the autosampler model | Agilent 1050 Agilent 1090 Agilent 1100 Agilent 1200 Agilent 7673A Agilent 7683A Agilent 79855(A) Agilent 5880 Agilent 5890 Agilent 6850. 27 Pos. tray Agilent 6890 Agilent CTC HTS/ HTC PAL Agilent CTC GC PAL Agilent 7693A AI 42 vial tray AI 60 vial tray AI |

