



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® Roll rim vial amber glass ND11, 2 ml,
12 x 32 mm, wide opening, with title block**

**€13.10
plus VAT &
Shipping**



Description

The neochrom rolling rim vials ND11, vol. 2 ml, 12 x 32 mm are made of borosilicate glass of the 1st hydrolytic class. The cooling point is 565 °C.

Suitable for:

- Agilent™ GC and LC autosamplers: EC-1000 to EC-1003
- Shimadzu autosamplers: EC-1000

Please order matching crimp caps separately.

Additional Information

No.	EC-1003
Manufacturer (Brand)	neochrom
EAN	4058072080266
Transport temperature	Room temperature
Volume	2 ml
Color bottle	brown
Height	32 mm
Floor outside	flat
Floor inside	flat
DM bottle	12 mm
DM Closure	11 mm
Neck size	far
Material bottle	Borosilicate glass (1st hydrolytic class)

