



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® ND20 headspace vials, 5 ml clear glass, 38 x20 mm, flat bottom , 100 pcs./pack

**€26.30
plus VAT &
Shipping**

Product Images



Description

The instruments on the market require a wide variety of bottles.

To help you identify the correct bottle, we have provided information on instrument suitability in our product overview. All bottles are made of 1st hydrolytic class glass

- A flat bottom may be necessary if the bottles need to run slightly downward in the instrument.
- Flat DIN rolled rim
A sealing disc has more contact surface on a flat DIN rolled rim.

Additional Information

No.	7-0798
Manufacturer (Brand)	neochrom
EAN	4058072030445
Transport temperature	Room temperature
Volume	5 ml
Color bottle	transparent (clear)
Height	38 mm
Floor outside	flat
Floor inside	flat
Neck size	far
DM bottle	20 mm
DM Closure	20 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
for the autosampler model	Varian 8400 100 Pos. Varian 8410-Autoinjector (10x2ml/ 6x5ml/ 5x10ml)

