



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® High Recovery vials with snap ring,
clear glass, wide opening, 2 ml, 12 x 32 mm;
Packing unit: 100 pcs**

**€96.60
plus VAT &
Shipping**



Description

neochrom High Recovery Vials with wide opening, Vol. 2ml, 12 x 32 mm guarantee a dead volume of less than 2 µl due to a specially shaped interior. This allows an autosampler to recover the entire sample. No glass inserts with springs are required. Matching lids must be ordered separately.

Additional Information

No.	EC-1281
Manufacturer (Brand)	neochrom
EAN	4058072079925
Transport temperature	Room temperature
Volume	2 ml
Color bottle	transparent (clear)
Height	32 mm
Floor outside	flat. with stand ring
Floor inside	flat
Neck size	far
DM bottle	12 mm
DM Neck	10 mm
DM Closure	12 mm
Material bottle	Borosilicate glass

