



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton snap lid vials w.lid vers.22 mm, 8ml PU 144 pcs. no.9-0222

€219.14
plus VAT &
Shipping

Product Images



Description

With convenient snap closure that seals tightly but opens easily.
Especially suitable for collecting and storing hygroscopic material.
Made of soda glass of hydrol. Class II.
Packed on trays.

Additional Information

No.	WD-0788
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0221
EAN	4058072427207
Transport temperature	Room temperature
Volume	8 ml
Color bottle	transparent (clear)
Height	39 mm
Closure / lid properties	closed
DM bottle	22 mm
Color closure/ lid	transparent (milky)
DM Closure	22 mm
Material bottle	Soda-lime glass (Wheaton 800. USP standard type 3)
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX bottle	510 C°
TBST MAX closure/ lid	80 °C
TBST MIN Closure/ lid	-100 °C

