



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wheaton snap lid vials w.lid vers.30mm, 16ml VE 144 pcs. no.9-0224

**€274.68**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

With convenient snap closure that seals tightly but opens easily.  
Especially suitable for collecting and storing hygroscopic material.  
Made of soda glass of hydrol. Class II.  
Packed on trays.

## Additional Information

---

No.	WD-0790
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0223
EAN	4058072427221
Transport temperature	Room temperature
Volume	16 ml
Color bottle	transparent (clear)
Height	40 mm
DM bottle	30 mm
DM Closure	30 mm
Color closure/ lid	transparent (milky)
Material bottle	Soda-lime glass (Wheaton 800. USP standard type 3)
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX bottle	510 °C
TBST MAX closure/ lid	80 °C
TBST MIN Closure/ lid	-100 °C
Closure / lid properties	closed

