



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton snap lid vials w.lid vers.30mm, 24ml, VE 144 pcs

€276.31
plus VAT &
Shipping

Product Images



Description

With convenient snap closure that seals tightly but opens easily.
Especially suitable for collecting and storing hygroscopic material.
Made of soda glass of hydrol. Class II.
Packed on trays.

Additional Information

No.	WD-0791
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0224
EAN	4058072427238
Transport temperature	Room temperature
Volume	24 ml
Color bottle	transparent (clear)
Height	52 mm
DM bottle	30 mm
DM Closure	30 mm
Color closure/ lid	transparent (milky)
Material bottle	Soda-lime glass (Wheaton 800. USP standard type 3)
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX bottle	510 °C
TBST MAX closure/ lid	80 °C
TBST MIN Closure/ lid	-100 °C
Closure / lid properties	closed

