



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® flat-bottom vials 1 ml clear glass,
35 x 7.8 mm, w. PE stopper 6 mm transp., 100
pcs./pack**

**€16.40
plus VAT &
Shipping**

Product Images



Description

A convenient to use and inexpensive bottle/closure combination for non-critical analyses made of glass of the 1st hydrolytic class.

A star-shaped diaphragm allows easy penetration of the PE stopper.

Recommended for use in the HPLC area.

For the 1 ml flat-bottomed tubes, stoppers with or without insertion lock for microinserts are available.

Stoppers without insertion barrier show a better valve effect with regard to the formation of a vacuum.

Additional Information

No.	7-0736
Manufacturer (Brand)	neochrom
EAN	4058072029906
Transport temperature	Room temperature
Volume	0.8 ml
Color bottle	transparent (clear)
Height	35 mm
Floor outside	flat
Floor inside	flat
Closure / lid properties	pierceable (with star) closed
DM bottle	7.8 mm
DM Closure	6 mm
Color closure/ lid	transparent (milky)
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	High-density polyethylene (HDPE)
for the autosampler model	Alcott 708 AL Alcott 728

