



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® flat-bottom vials 1 ml, clear glass  
with 8 mm PE stopper transparent, 40x8.2 mm,  
100 pcs./pack**

**€13.40  
plus VAT &  
Shipping**



## Description

A convenient to use and inexpensive bottle/closure combination for non-critical analyses made of glass of the 1st hydrolytic class.

A star-shaped diaphragm allows easy penetration of the PE stopper.

Recommended for use in the HPLC area.

For the 1 ml flat-bottomed tubes, stoppers with or without insertion lock for microinserts are available.

Stoppers without insertion barrier show a better valve effect with regard to the formation of a vacuum.

## Additional Information

No.	7-0737
Manufacturer (Brand)	neochrom
EAN	4058072029913
Transport temperature	Room temperature
Volume	1 ml
Color bottle	transparent (clear)
Height	40 mm
Floor outside	flat
Floor inside	flat
Closure / lid properties	pierceable (with star) closed
DM bottle	8.2 mm
Color closure/ lid	transparent (milky)
DM Closure	8 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	High-density polyethylene (HDPE)
for the autosampler model	Agilent 7693A Cecil Instruments AutoQuest Dionex AS 40 Shimadzu LC-2010A/C 350 pos. tray Shimadzu SIL-10A/ AF/ AP/ Ai/ AxL Rack MTP2 192 pos. Shimadzu SIL-10ADvp Shimadzu SIL-10HTA/ HTC 350 pos. tray Shimadzu SIL-20A/ AC (Prominence) 175 vial tray Waters®

