



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® flat-bottom vials 2 ml clear glass,
31.5 x 11.6 mm, PE stopper, 100 pcs./pack**

**€14.00
plus VAT &
Shipping**

Product Images



Description

A convenient to use and inexpensive bottle/closure combination for non-critical analyses made of glass of the 1st hydrolytic class.

A star-shaped diaphragm allows easy penetration of the PE stopper.

Recommended for use in the HPLC area.

For the 1 ml flat-bottomed tubes, stoppers with or without insertion lock for microinserts are available.

Stoppers without insertion barrier show a better valve effect with regard to the formation of a vacuum.

Additional Information

No.	7-0742
Manufacturer (Brand)	neochrom
EAN	4058072029951
Transport temperature	Room temperature
Volume	1.5 ml
Color bottle	transparent (clear)
Height	31.5 mm
Floor outside	flat
Floor inside	flat
Closure / lid properties	pierceable (with star) closed with fixation possibility for micro inserts
DM bottle	11.6 mm
Color closure/ lid	transparent (milky)
DM Closure	12 mm
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	High-density polyethylene (HDPE)
for the autosampler model	Alcott 738 Carlo Erba AS100 Carlo Erba A200LC Carlo Erba AS300 Fisons AS100 Fisons A200LC Fisons AS300 GBC Avanta Ultra Z Gilson 201/202 Gilson 221/222 Gilson 231/401 Gilson 232/402 Gilson Aspec Shimadzu LC-2010A/C 140 Pos. Tray Shimadzu SIL-6B/ 7A/ 8A/ 9

