



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® flat-bottomed vials 2 ml amber
glass, 31.5 x 11.6 mm, PE stopper, 100 pcs./pack**

**€15.20
plus VAT &
Shipping**

Product Images



Description

A convenient to use and inexpensive bottle/closure combination for non-critical analyses made of glass of the 1st hydrolytic class.

A star-shaped diaphragm allows easy penetration of the PE stopper.

Recommended for use in the HPLC area.

For the 1 ml flat-bottomed tubes, stoppers with or without insertion lock for microinserts are available.

Stoppers without insertion barrier show a better valve effect with regard to the formation of a vacuum.

Additional Information

| | |
|---------------------------|---|
| No. | 7-0743 |
| Manufacturer (Brand) | neochrom |
| EAN | 4058072029968 |
| Transport temperature | Room temperature |
| Volume | 1.5 ml |
| Color bottle | brown |
| Height | 31.5 mm |
| Floor outside | flat |
| Floor inside | flat |
| Closure / lid properties | pierceable (with star) closed with fixation possibility for micro inserts |
| DM bottle | 11.6 mm |
| DM Closure | 12 mm |
| Color closure/ lid | transparent (milky) |
| Material bottle | Borosilicate glass (1st hydrolytic class) |
| Material closure/ lid | High-density polyethylene (HDPE) |
| for the autosampler model | Alcott 738 Carlo Erba AS100 Carlo Erba A200LC Carlo Erba AS300 Fisons AS100 Fisons A200LC Fisons AS300 GBC Avanta Ultra Z Gilson 201/202 Gilson 221/222 Gilson 231/401 Gilson 232/402 Gilson Aspec Shimadzu LC-2010A/C 140 Pos. Tray Shimadzu SIL-6B/ 7A/ 8A/ 9 |

