



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657




---

**Lenz Chromatography Column with grooves  
and PTFE stopcock, eff. length 200 mm, 35 ml,  
flared rim**

**€54.68  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Simple chromatography columns with punctures above the stopcock for placing a cotton swab with flared rim with flared rim with PTFE plug and stopcock safety device; made of DURAN® tube.

## Additional Information

---

No.	LZ-1398
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3401
EAN	4058072311377
Transport temperature	Room temperature
diameter	15 mm
Material	Borosilicate glass
Volume	35 ml
Length	200 mm
Standard ground sleeve	14/23
Type cock / valve	PTFE - tap

