



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz chromatography column with punctures
and PTFE stopcock, eff. length 400 mm, 125 ml,
flared rim**

**€60.07
plus VAT &
Shipping**

Product Images



Description

Simple chromatography columns with punctures above the stopcock for placing a cotton swab with flared rim with flared rim with PTFE plug and stopcock safety device; made of DURAN® tube.

Additional Information

No.	LZ-1399
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3402
EAN	4058072311384
Transport temperature	Room temperature
diameter	20 mm
Material	Borosilicate glass
Volume	125 ml
Length	400 mm
Standard ground sleeve	14/23
Type cock / valve	PTFE - tap

